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Version 01 – 17.Mai 2021

Meoplast® Medical Implant System

Reprocessing instructions for the Meoplast Medical Implant System

DISCLAIMER

The products mentioned here are components of a comprehensive treatment concept and may only be used in combination with the associated original products in accordance with the instructions and recommendations of Meoplast® Medical GmbH. The non-recommended use of third-party products in combination with Meoplast products will void the warranty and any other express or implied obligations of Meoplast will be void. The user of Meoplast products must determine whether the product is suitable for a particular patient under the given conditions. Meoplast assumes no liability, express or implied, for direct or indirect damages, penalties including damages or other damages arising from or in connection with errors in professional judgment or practice in the context of the use of Meoplast products. The user is also obliged to regularly inform himself about the latest developments with regard to Meoplast products and their application. If in doubt, contact Meoplast. Since the use of the products is under the control of the user, the user assumes responsibility. Meoplast accepts no liability whatsoever for damage resulting from the use of the product.

SCOPE OF APPLICATION OF THE INSTRUCTIONS FOR USE

These instructions for use are valid for all components of the Meoplast implant system that may be reprocessed (reusable instruments) in accordance with their intended purpose and labeling, taking into account the Medical Devices Regulation 2017/745. Reusable products include all insertion tools of the Meoplast Medical Implant System (drill extension, hexagonal instruments, implant insertion instruments, torque ratchet, depth stops and surgical trays). The implants, abutments and impression tools are for single use and are not designed for reuse.

SAFETY

For reasons of hygiene and health protection, all reusable medical devices distributed by Meoplast Medical (see scope of the instructions for use) must be cleaned, disinfected and sterilized before/after each application to the patient or in the patient's mouth. This

also applies to the first use after receipt of the products, which are delivered non-sterile and must be sterilized.

The hygiene regulations of the individual national and international legal provisions of the dentistry/medical practices, the hospitals as well as the dental laboratories must be complied with.

For cleaning, disinfection and sterilization, only devices that can demonstrate a validated procedure may be used. The following meters for cleaning and sterilization must be complied with, as they have been subjected to a validated process. When carrying out the preparation with different parameters, the cleanliness of the product cannot be guaranteed.

All Meoplast Medical implants are intended for single use and are delivered sterile. These must not be cleaned, disinfected and sterilized again: they are marked as follows:

Do not re-sterilize

Do not reuse

Sterilized by irradiation

INSTRUCTIONS FOR CLEANING AND STERILIZING

Any reprocessing must be carried out by trained personnel in a room specially equipped for this purpose (with a clean and unclean zone).

Successful sterilization is preceded by efficient cleaning and disinfection. Before/after each application to the patient, the instruments and prosthetic components must be cleaned, disinfected and sterilized.

Dirty/used instruments must be separated from the (still) clean instruments during use. These are not to be put back into the surgical tray. After successful cleaning and disinfection, the instruments are sorted back into the surgical tray and then sterilized as a set.

Manual pre-cleaning

5. Immediately after use, place all components to be cleaned in a liquid instrument cleaner for 30

minutes (e.g. INSTRU PLUS liquid cleaner for instruments and endoscopes from Dr. Schumacher, concentration 5mg/L). If the components can be dismantled, they must first be disassembled into their individual parts.

6. Remove any other visible residue with a nylon instrument brush. Cavities, bushings and constrictions must be treated twice.
7. Rinse all components under cold, running water (drinking water quality).
8. Visual inspection of the purity of all individual parts. Magnifying glasses can be added as a support.

Repeat steps 1 through 4 if residue is still visible on the components.

Machine cleaning and disinfection

5. Place all components in an appropriately approved (small) parts holder in the washer-disinfector (according to ISO 15883). Previously dismantled components must be reassembled.

All products must be arranged in such a way that they can be hit directly by the spray.

6. Cleaning and disinfection must be started using the Vario TD program and a cleaning agent for mechanical cleaning (e.g. Neodisher MediClean Dental from Dr. Weigert). Other cleaning programs are not covered by the validation.
7. Drying of the instruments with purified compressed air. Pay attention to hard-to-reach areas such as bushings.
8. Visual inspection of the purity and integrity of all components. Magnifying glasses can be added as a support.

Repeat steps 1 through 4 if residue is still visible on the components.

Packaging and sterilization

All components supplied non-sterile must not be sterilized in their original packaging.

5. Assemble the surgical tray with the appropriate components. All instruments are to be classified accordingly. V Then open up the tray.
6. The surgical tray (if applicable, the individual components that are not to be sorted into the surgical tray) is packed in a sterile bag in accordance with ISO 11607-1 and sealed by means of a sealing machine.

Care must be taken to ensure that the products in the sterile bag are not energized. This could, among other things, lead to damage to the sterile barrier and thus to the loss of sterility.

7. Sterilization is carried out with a class B small steam sterilizer with fractional backing vacuum in accordance with ISO 13060. Sterilization parameters other than those mentioned below are not covered by the validation. Thus, no guarantee can be given that the "sterile" status will be achieved.

Temperature	Pressure	Holding time	Dry season
134°C ± 1%	3 bar ± 5%	≥ 3 min.	≥ 5 min.

8. The marking of the reprocessed, packaged components must contain the following information:
 - a. Name of the product(s);
 - b. Sizes
 - c. Release details
 - d. Release decision
 - e. Sterilization cycle and sterilization date
 - f. Expiry date and sterile goods storage period

The sterilized components should not be stored for more than 6 days. It is recommended to carry out sterilization before application to the patient.

DISPOSAL

Products that have a loss of function and/or corrosion points should be removed from the application immediately and no longer used on the patient.

Defective products must be disposed of in accordance with local and environmental regulations, taking into account the level of contamination in question.

Manufacturer:

Meoplant Medical GmbH
 Malchiner Straße 99 – 12359 Berlin / Germany
 Tel. +49 30 809334166
 E-mail: info@meoplant.de
 Web: https://www.meoplant.de

LEGEND

Attention: Pay attention to accompanying documents

Observe the instructions for use

Batch number (see packaging)

Article (see packaging)

CE marking for Class I products

0297
 CE marking for class IIa/IIb products

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–All rights reserved–